



FIRSTSAVE® AED G3

automated external defibrillator



FEATURES

- Simple, easy to use, one-button operation
- Patented technology performs daily, weekly, and monthly self-tests
- More instructive voice prompts guide user through rescue
- Outstanding voice prompt quality and clarity
- Lightweight in a new user-friendly form factor

TECHNOLOGY AND PERFORMANCE — IN A VALUE PRICED DESIGN

Introducing the Cardiac Science FirstSave AED G3.

The next generation FirstSave AED has arrived. Technology and value – the FirstSave AED G3 offers customers our leading edge AED technology in a value priced AED design.

The new FirstSave AED G3 features our daily, weekly, and monthly self-testing all automatically testing three critical components (pads presence, software/electronics, and battery system). The FirstSave also incorporates our advanced technology features such as the Company's patented and programmable RHYTHMx® analysis software and STAR® biphasic defibrillation energy waveform.



One-button operation makes the FirstSave AED G3 easy to use and simple to operate.



To learn more about the best choice in AEDs, visit our website www.cardiacscience.com, or email info@cardiacscience.com, or call 1.866.289.5649 (USA and Canada only).



CARDIAC SCIENCE

SPECIFICATIONS

DEFIBRILLATOR

Operation	Semi-automatic
Waveform	Biphasic truncated exponential
Energy (J) range available	200VE, 200VE, 200VE*
Shock time cycle (Typical)	Less than 55 seconds for three shock series
Controls	One-button operation
Voice prompts	Comprehensive voice instructions guide user through rescue process
Visible indicators	SmartGauge Battery Status Indicator, Service Indicator, Pads Indicator Audible alerts Voice Prompt, System Alert Built in automatic synchronization
Synchronized shock	Yes
Pacemaker pulse detection	No
Programmable	Yes
Pediatric capability	Yes

PADS

Minimum combined surface area	228cm ²
Extended length of lead wire	1.3m
Supplied	Pre-connected to the AED, test for presence
Type	Adult, pre-gelled, self-adhesive, disposable, non-polarized (identical pads can be placed in either position) defibrillation pads
Shelf life	2 years

*USA only.

BATTERY

Type	IntelliSense [®] lithium battery
Warranty	1 year, full operational replacement

AUTOMATED SELF-TESTS

Daily	Battery, pad presence, internal electronics, no energy charge, SHOCK/CONTINUE button, and software
Weekly	Battery, pad presence, internal electronics, partial energy charge, SHOCK/CONTINUE button, and software
Monthly	Battery, pad presence, internal electronics, full energy charge cycle, SHOCK/CONTINUE button, and software

EVENT DOCUMENTATION

Type	Internal memory
Internal memory	34 minutes ECG data with event annotation
Playback	Viewable via RescueLink [®] software via PC
Communications	Serial port or USB (via adapter) for PC with Windows
Clock synchronization	Precise synchronization of AED time to dispatch time

PHYSICAL DIMENSIONS

Height	8 cm (3.3 in)
Width	27 cm (10.6 in)
Depth	31 cm (12.4 in)
Weight	3.10 kg (6.6 lb)

ORDER INFORMATION

MODEL

9300C-201

DESCRIPTION

FirstSave AED G3 Automated External Defibrillator

Each AED package includes: (1) defibrillator, (1) IntelliSense battery (9143), (1) pair of defibrillation pads, and (1) Quick Start Tool Kit: includes Quick Start Guide, CD-ROM with AED Manual, Training Video, RescueLink and MDLink, and serial communication cable.